



# State Police Face Mask

www.tuttomilitare.com



www.tuttomilitare.com

www.tuttomilitare.com

www.tuttomilitare.com

**Product price:**

**€9.50 tax included**

**Product codes:**

Reference: Art. 2969

EAN13: 0658606579460

UPC: 658606579460

**Product features:**

## Product short description:

Official Polizia di Stato face mask for those who want a practical, understated accessory that conveys authority and respect. Made from high-quality materials, it offers comfortable wear and reliable protection for everyday use.

- **Material:** High-quality fabric for comfort and durability.
- **Design:** Elegant, distinctive look with the Polizia di Stato logo.
- **Use:** Ideal for daily wear, commuting and public settings.
- **Fit:** Comfortable, secure fit with adjustable elastic ear loops.
- **Symbol:** A sign of commitment and dedication.

## Product description:

## Product overview

The Mascherina Polizia di Stato is a duty face mask designed for Polizia di Stato personnel. It combines a washable, elastic technical fabric with the official Polizia di Stato logo to provide comfortable daily protection for police duty and tactical use.



Definition: A reusable, washable protective mask for law enforcement use, produced to meet the requirements cited in Decreto Legge n.18 (17/03/2020) and the related Ministry of Health circular.

### Main features

Built from a certified elastic, non-run fabric, this mask prioritizes fit, breathability and clear identification for duty use.

- **Official logo:** Printed with the golden eagle and the RI monogram, labeled "Polizia di Stato".
- **Material:** Elastic, non-run fabric — 78% PES and 22% EA, Oeko-Tex Standard 100 certified.
- **Washable and reusable:** Durable through up to 20 washes; retains hydrophobic treatment and breathability when washed as recommended.
- **Comfort and fit:** High elastic modulus and ergonomic shape ensure close fit for extended wear on duty.
- **Breathability:** Fabric is breathable and pleasant to the touch for long shifts.
- **Laboratory observation:** SEM analysis from the University of Padova found fabric pore dimensions on the order of 50 microns, indicating potential to reduce droplet spread.

### How to use / When it's useful

Use the Mascherina Polizia di Stato as a duty mask during public interactions, patrols and routine operations where a washable, identifiable face covering is required. It is intended for personnel who must display official identification on their protective equipment.

Quick answer: Wear during shifts when a face covering is required by service rules. In most cases, follow your unit's PPE guidance and laundering procedures to keep the mask effective.

### Care and notes

Simple care keeps performance and appearance intact; follow laundering instructions and medical or departmental guidance.

- Wash up to 20 times at 90°C as specified; this preserves hydrophobic treatment and breathability.
- Do not alter the mask or remove official insignia.
- Sale is restricted to Polizia di Stato personnel; purchase requires ID verification (badge copy or service attestation) or registering the Questura delivery address for verification.
- The mask is intended to help contain droplet spread but is not labeled here as a medical-grade respirator; follow operational PPE protocols.

### FAQ

**Q:** Who can buy the Mascherina Polizia di Stato?



**A:** Purchase is reserved for Polizia di Stato personnel. Verification requires a copy of the service badge or an official service attestation, or a verified Questura delivery address.

**Q:** How many times can the mask be washed?

**A:** The mask can be washed up to 20 times at 90°C while maintaining its hydrophobic and breathable properties when laundered as recommended.

**Q:** What is the fabric composition?

**A:** The mask is made from elastic, non-run fabric: 78% PES (polyester) and 22% EA (elastane), certified Oeko-Tex Standard 100.

**Q:** Does the mask meet current legal guidance?

**A:** The product is produced in line with the indications of Decreto Legge n.18 (17/03/2020) and the related Ministry of Health circular referenced in the source text.

**Q:** Does the fabric stop infectious droplets?

**A:** SEM analysis by the University of Padova showed fabric porosities around 50 microns; this suggests the material can help reduce droplet spread, but it should be used alongside other infection-control measures.

**Q:** How is identity verified for purchase?

**A:** Verification is done by submitting a badge copy or service attestation via the seller's contact process, or by listing the Questura as the delivery address during registration so the address can be checked before allowing purchase.

**Tip:** Choose the correct variant and ensure your delivery details match your Questura to speed verification. Browse related duty and tactical gear to complete your issued equipment while following your unit's PPE guidance.

**Product attributes:**

**Product gallery:**

