



German Red Cross 2nd Class Ribbon



Product price:

€3.50 tax included

Product codes:

Reference: Art. 3011

EAN13: 0658606579880

UPC: 658606579880

Product features:

Product short description:

The Nastrino Croce Rossa Tedesca II Classe is a ribbon honoring humanitarian and medical service. Suitable for collectors, reenactors and anyone displaying or wearing period medical insignia.

- **Material:** Quality materials with an elegant, understated design for display or wear on uniforms and formal attire.
- **Purpose:** Recognition of service in medical and humanitarian support; a symbol of dedication and compassion.
- **Wear:** Fits on uniforms or formal dress; also suitable for display in collections.
- **Ideal for:** Collectors, history enthusiasts and institutions preserving Red Cross memorabilia.

Product description:

Product overview

The Nastrino Croce Rossa Tedesca II Classe is a service ribbon issued to civilian or military personnel who have rendered honorable services to the Croce Rossa Tedesca. It is intended for uniform wear, collectors, military surplus sellers, and display with tactical or duty/police kit.



Definition: a ribbon bar representing the Croce Rossa Tedesca Second Class award, suitable for mounting on ribbon bars or purchasing as loose ribbon for medal work.

Main features

Clear, regulation-style ribbon available in multiple formats and attachment options for collectors and uniform restorers.

- **Accurate design:** 37 mm wide, white field with an 18 mm central vertical red stripe and two 2 mm silver stripes at the edges.
- **Materials:** Made from durable, high-quality ribbon fabric for long-lasting appearance.
- **Formats:** Regulation flat version or slightly domed Elite variant.
- **Attachment options:** Rear fittings available as clamps or screw-type fasteners (select in the dropdown).
- **Loose ribbon available:** Sold by the metre with a minimum purchase of 50 cm for custom medal work.
- **Suitable for ribbon bars:** Can be mounted on multi-place ribbon bars (two or more positions).

How to use / When it's useful

Quick answer: Wear on the uniform where service ribbons are displayed or use the loose ribbon for medal mounting and restoration work. Select clamp or screw fittings based on your ribbon bar or uniform requirements.

In most cases, collectors and military surplus retailers buy the regulation flat form, while restorers and display builders may prefer the Elite domed version or loose ribbon to cut and sew. Choose the attachment type that matches your ribbon bar hardware.

Care and notes

Basic care will keep the ribbon looking correct for display and wear.

- Avoid harsh cleaning chemicals; spot clean with mild soap and water if needed.
- Store flat or on a ribbon bar to prevent creasing, especially for the Elite domed variant.
- Minimum order for loose medal ribbon is 50 cm—select the appropriate option in the dropdown when ordering.
- Check your selected rear attachment (clamp or screw) before finalizing the order to ensure compatibility with your ribbon bar.

FAQ

Q: Who can receive the Nastrino Croce Rossa Tedesca II Classe?



A: It is awarded to civilian or military personnel who have rendered honorable services to the Croce Rossa Tedesca. This listing offers the ribbon for collectors and uniform use.

Q: What are the ribbon dimensions?

A: The ribbon is 37 mm wide, with an 18 mm central vertical red stripe and two 2 mm silver edge stripes, matching the described pattern.

Q: Which formats are available?

A: Available as the regulation flat ribbon or the slightly domed Elite variant. You can also buy the ribbon loose by the metre (minimum 50 cm).

Q: What attachment options are offered?

A: Rear fittings are available as clamps or screw-type fasteners. Please select the desired attachment from the product dropdown when ordering.

Q: Can this ribbon be mounted on a multi-place ribbon bar?

A: Yes. The ribbon can be mounted on two-or-more place ribbon bars designed to hold multiple ribbons on a single metal base.

Q: Is the ribbon always in stock?

A: Demand can be high and availability may be limited. If you need a specific format or attachment, order promptly to secure the correct variant.

Tip: Choose the correct format and rear attachment for your ribbon bar or display, and consider the loose ribbon option (minimum 50 cm) for custom medal work. Browse related ribbon bars and mounting options to complete your display or uniform setup.

Product attributes:

Nastrini : Formato 37 mm., Confezione Regolare, Confezione Elite, Confezione Regolare con Attacco a Clip, Confezione Regolare con Attacco a Vite, Confezione Elite con Attacco a Clip, Confezione Elite con Attacco a Vite, 50 cm. Nastro 37 mm.

Product gallery:

