



# Second Lieutenant Medical Officer Cap

www.tuttomilitare.com



www.tuttomilitare.com

www.tuttomilitare.com

www.tuttomilitare.com

## Product price:

**€95.00 tax included**

## Product codes:

Reference: Art. 3043

EAN13: 0658606580206

UPC: 658606580206

## Product features:

## Product short description:

The Berretto Sottotenente Medico is a refined officer's cap for medical personnel and anyone seeking a polished, professional accessory. Made from high-quality materials, it delivers comfort and lasting wear while conveying prestige and competence.

- **Material:** High-quality construction for comfort and durability.
- **Design:** Elegant, distinctive look with refined details that signal professionalism.
- **Use:** Ideal for medical officers, formal occasions, or adding a smart touch to everyday wear.
- **Care:** Easy to maintain; follow standard garment care to preserve appearance.

## Product description:

## Product overview

This is the Berretto Sottotenente Medico: a formal service cap designed for medical officers serving in the Esercito. It is intended for active military personnel who require the correct rank insignia and a durable, service-ready cap.



**Definition:** A cordellino kaki cap with damask band and gold-wired embroidered badge for the Sottotenente Medico rank. The cap combines traditional appearance with practical materials used in tactical and duty uniforms.

### Main features

Key details and materials that define the cap.

- **Material:** Cordellino kaki fabric with a damask (damascata) band.
- **Embroidered badge:** Gold canutiglia embroidery on kaki cloth showing a five-point star topped by a turreted crown with two caducei and a central white disc with a red cross.
- **Rank fitting:** Supplied with the Soggolo for the Sottotenente rank; Tenente and Capitano variants available via the dropdown menu.
- **Packaging time:** Cap is packaged within 3 days from order.
- **Intended use:** Reserved for active military personnel; sale requires identity/command verification.

### How to use / When it's useful

Use this cap as part of formal uniform duties, ceremonies, or when a compliant rank cap is required. It suits medical officers who need a duty-cap that matches Esercito regulations and appearance standards.

**Quick answer:** Use it for official duty, parades, and formal occasions where the Medico cap is required. In most cases, select the correct rank variant (Sottotenente, Tenente, Capitano) from the product options before ordering.

### Care and notes

Practical care tips and important purchase notes.

- Store in a dry place to preserve shape and embroidery; avoid crushing the cap.
- Spot clean or consult a professional cleaner for embroidered and gold-wired details.
- Sale reserved to active military personnel: purchase requires a copy of your service ID or delivery to your Command headquarters for verification.
- Limited availability: items for military supply and military surplus or tactical duty use are subject to verification before purchase is approved.

### FAQ

**Q:** Who can buy the Berretto Sottotenente Medico?

**A:** Only active military personnel may purchase this cap. Verification is required by providing a copy of your service ID or indicating your Command headquarters as the delivery address.



**Q:** How long does it take to prepare the cap for shipment?

**A:** The cap is packaged within 3 days from the order date. Verification steps may affect the final shipping time.

**Q:** Can I order the cap in other ranks?

**A:** Yes. The cap is available in Sottotenente, Tenente, and Capitano variants—select the rank from the dropdown menu when ordering.

**Q:** What exactly is the badge design?

**A:** The badge is a five-point star surmounted by a turreted crown with two caducei; a white disc with a red cross is centered on the badge. The embroidery uses gold canutiglia on kaki cloth.

**Q:** What materials are used for the cap?

**A:** The cap is made in cordellino kaki fabric with a damask band and a gold-wired embroidered badge on kaki cloth, designed for durability in duty and tactical use.

**Q:** How is membership in the Corps verified?

**A:** Verification is done by sending a copy of your service ID or by specifying your Command headquarters as the delivery point; the Command details will be checked before purchase approval.

**Tip:** Choose the correct rank variant before adding the cap to your cart and check the packaging time. Consider browsing related uniform and tactical duty items to complete your set while respecting the verification requirements.

**Product attributes:**

Misure per Berretti e Baschi: 54 , 55, 56, 57, 58, 59, 60

Berretti con Soggoli Ufficiali: Sottotenente, Tenente , Capitano, Maggiore, Tenente Colonnello, Colonnello

**Product gallery:**



www.tuttomilitare.com



www.tuttomilitare.com

www.tuttomilitare.com

www.tuttomilitare.com

