



Italian Red Cross Time of Kindness Excellence Ribbon



Product price:

€5.50 tax included

Product codes:

Reference: Art. 3076

EAN13: 0658606580534

UPC: 658606580534

Product features:

Product short description:

The Nastrino di Eccellenza Tempo della Gentilezza Croce Rossa Italiana is a ribbon badge that honors exceptional humanitarian service. It is intended for volunteers, staff and supporters of the Italian Red Cross who wish to carry a visible sign of solidarity and compassion.

- **Material:** Crafted from quality materials for a respectful, durable finish.
- **Use:** Recognition ribbon suitable for wear on uniforms or for display.
- **Meaning:** A tangible symbol of dedication, altruism and significant community contribution.
- **Ideal for:** Volunteers, Croce Rossa Italiana personnel, donors and those commemorating acts of kindness.

Product description:

Product overview

The Nastrino di Eccellenza "Il Tempo della Gentilezza" Croce Rossa Italiana is a service ribbon designed for Croce Rossa Italiana (CRI) volunteers and staff who served during the Covid-19 emergency. It is aimed at CRI personnel, units assembling uniform decorations, and collectors in tactical, military surplus and duty/police memorabilia channels.



Definition: this is a service ribbon (nastrino) awarded to CRI personnel for service, coordination or notable effort during the Covid-19 health emergency.

Main features

Clear visual design, durable materials and multiple mounting options make this ribbon suitable for uniform wear or display.

- **Distinctive design:** white and red Croce Rossa ribbon crossed at the center by four alternating yellow and black staves.
- **Five-point silver star:** a small silver star mounted at the ribbon center for added prestige.
- **High-quality materials:** made to resist wear and preserve appearance over time.
- **Available formats:** regulation flat format and slightly domed Elite variant.
- **Attachment options:** single rear fitting with clamps or screws; selectable from the product options.
- **Mounting choices:** offered as a loose ribbon for medal use or fitted on multi-ribbon bars (two or more positions) on a single metal base.

How to use / When it's useful

Quick answer: wear the ribbon on the uniform ribbon bar or display it on a medal or multi-ribbon holder according to CRI practice. In most cases, choose the regulation flat format for standard uniform wear or the Elite variant for a slightly raised look.

Select the rear attachment (clamps or screws) when ordering to match your mounting needs. For display or collections, order the ribbon as loose for medal mounting or add it to a multi-ribbon base for organized presentation.

Care and notes

Basic care preserves the ribbon's appearance and symbolic value; choose the correct variant and attachment before purchase.

- Store in a dry place and avoid prolonged exposure to direct sunlight to limit fading.
- Clean gently with a soft, dry cloth; avoid solvents and strong detergents that can damage the finish.
- Available in limited quantities—check variant and attachment options before completing your order.
- Select the attachment type in the product options to ensure correct mounting on your ribbon bar or medal.

FAQ

Q: Who can receive the Nastrino di Eccellenza "Il Tempo della Gentilezza"?



A: The ribbon is awarded to Croce Rossa Italiana volunteers and employees who served during the Covid-19 emergency, including those who provided coordination or notable service.

Q: What formats and sizes are available?

A: The ribbon is sold in the regulation flat format and in a slightly domed Elite variant. It is also available as a loose ribbon for medal mounting.

Q: How can I personalize or mount the ribbon?

A: Choose rear attachments (single fitting with clamps or screws) from the product options. You can also order the ribbon mounted on multi-ribbon bars with two or more positions.

Q: What do the yellow and black staves mean?

A: The yellow/black staves reference maritime signal colors historically used to indicate quarantine or special caution; here they underline the ribbon's link to epidemic response.

Q: How do I buy this ribbon?

A: Add the preferred format and attachment to your cart and complete the checkout online. Check available stock and chosen options before finalizing the purchase.

Q: Is this an official CRI award and how can I confirm eligibility?

A: The product represents the Nastrino di Eccellenza as described for CRI personnel who served during Covid-19. For formal confirmation of conferment or eligibility, please contact Croce Rossa Italiana or your unit administration.

Tip: Choose the correct variant (flat or Elite) and the appropriate rear attachment before ordering, and consider a multi-ribbon holder if you plan to display multiple service ribbons; browse related ribbons and mounting options to complete your set.

Product attributes:

Nastrini : Formato 37 mm., Confezione Regolare, Confezione Elite, Confezione Regolare con Attacco a Clip, Confezione Regolare con Attacco a Vite, Confezione Elite con Attacco a Clip, Confezione Elite con Attacco a Vite, 50 cm. Nastro 37 mm., 20 cm. di Nastro 37 mm.

Product gallery:

