




# 2002 Campobasso-Foggia Earthquake Commemorative Ribbon Red Cross



**Product price:**

**€3.50 tax included**

**Product codes:**

Reference: -

EAN13: 0658606553880

UPC: 658606553880

**Product features:**

## **Product short description:**

The "Nastrino Terremoto Campobasso Foggia 2002" is a commemorative ribbon honoring the Italian Red Cross response to the 2002 earthquake. It's aimed at collectors, supporters, and anyone who wants to pay tribute to the volunteers who assisted the affected communities.

- **Material:** Produced with high-quality materials as a tangible symbol of solidarity.
- **Use:** Wear or display to remember the relief efforts and show support for disaster response.
- **Tribute:** Commemorates the intervention of the Croce Rossa Italiana in Campobasso and Foggia, 2002.
- **Audience:** Suitable for collectors, humanitarian supporters, and exhibition displays.

## **Product description:**

## **Product overview**

The Nastrino Terremoto Campobasso Foggia 2002 Croce Rossa is a commemorative ribbon for collectors, uniform wearers, Red Cross volunteers and anyone honoring relief and reconstruction work after the 2002 earthquakes in Campobasso and Foggia. It is suitable for military surplus collectors, tactical and police uniform displays, and memorial presentations.



Definition: this ribbon was established by Italian Red Cross (CRI) ordinance on 9 December 2005, n. 1599 as a symbol of solidarity and reconstruction for communities affected by seismic events.

## Main features

Key features and options at a glance.

- **Standard size:** 37 mm width, matching common ribbon bar standards.
- **Two finishes:** flat regulatory format or slightly domed "Elite" variant.
- **Mounting options:** adjustable rear attachment with clamps or screws.
- **Loose ribbon option:** available as a separate ribbon strip for medal mounting.
- **Can be mounted:** suitable for single holders or multi-ribbon bars (porta nastrini).
- **Official origin:** instituted by CRI ordinance (9 December 2005, n. 1599).

## How to use / When it's useful

Quick answer: wear the ribbon on a ribbon bar, add it to a medal ribbon, or use it in a display to commemorate aid and reconstruction efforts. It fits standard metal ribbon mounts and multi-place holders.

In most cases, choose the flat format for a traditional look or the Elite domed variant for a slightly raised appearance. Select clamps or screws for the rear attachment based on the garment or display plate thickness to ensure a secure fit.

## Care and notes

Basic care and selection notes to keep the ribbon presentable and correctly mounted.

- Store dry and avoid prolonged exposure to direct sunlight to preserve color and finish.
- Check the chosen rear attachment (clamp or screw) before wearing to prevent loss or damage.
- Size is 37 mm; confirm compatibility with your ribbon bar or medal before purchase.
- Availability is limited and demand can be high; consider your preferred format early.

## FAQ

**Q:** What are the dimensions of the Nastrino Terremoto Campobasso Foggia 2002?

**A:** The ribbon follows the standard 37 mm width. This size matches common ribbon bars and mounting plates.

**Q:** What mounting options are available?



**A:** You can choose a rear attachment with clamps or screws. Pick the option that best suits your garment or display plate thickness.

**Q:** Are there different formats of the ribbon?

**A:** Yes. It is offered in the flat regulatory format and a slightly domed "Elite" variant for a raised look.

**Q:** Can I place this ribbon on a multi-ribbon holder?

**A:** Yes. The ribbon is compatible with two-place or multi-place metal ribbon holders (porta nastri) when mounted correctly.

**Q:** Who issued this ribbon and when was it established?

**A:** The ribbon was instituted by the Italian Red Cross (CRI) by ordinance on 9 December 2005, n. 1599 to commemorate the 2002 earthquakes in Campobasso and Foggia.

**Q:** Is there a loose ribbon option for medals?

**A:** Yes. A loose ribbon strip is available for those who wish to attach the ribbon to existing medals or custom mounts.

**Tip:** Choose the format (flat or Elite) and attachment (clamp or screw) that match your uniform or display. Check the 37 mm size against your ribbon bars and browse related commemorative ribbons and holders to complete your display.

#### **Product attributes:**

Nastri : Formato 37 mm., Confezione Regolare, Confezione Elite, Confezione Regolare con Attacco a Clip, Confezione Regolare con Attacco a Vite, Confezione Elite con Attacco a Clip, Confezione Elite con Attacco a Vite, 50 cm. Nastro 37 mm.

#### **Product gallery:**

