



CRI Volunteer Nurses Ribbon



Product price:

€3.50 tax included

Product codes:

Reference: -

EAN13: 0658606553927

UPC: 658606553927

Product features:

Product short description:

Nastrino Infermiere Volontarie CRI is the service ribbon for volunteer nurses of the Croce Rossa Italiana. It is designed to be worn with pride on uniforms and official attire.

- **Material:** Made from high-quality materials with neat stitching for a durable, professional finish.
- **Use:** Intended for display on uniforms and official clothing to identify volunteer nurse affiliation.
- **Design:** Elegant and distinctive, it pays homage to the tradition and dedication of CRI volunteer nurses.
- **Purpose:** A meaningful accessory to express support, commitment and belonging to a humanitarian cause.

Product description:

Product overview

The **nastrino Infermiere Volontarie CRI** is a service ribbon intended for volunteer nurses of the Croce Rossa Italiana, used to display affiliation and service on uniforms or display mounts. One simple definition: it is a 37 mm regulation ribbon bar representing the role of a CRI volunteer nurse.



Available in a flat regulation finish or a slightly domed "Elite" variant, the ribbon can be ordered with an optional ribbon holder that attaches by rear clamp or by screws. This item suits uniform, duty, tactical, and military surplus collections or displays.

Main features

Compact, regulation-compliant, and selectable with different mounting options to match daily use or collection display needs.

- **Regulation size:** 37 mm standard format.
- **Two finishes:** flat regulation version and slightly domed "Elite" variant.
- **Optional holder:** choose rear clamp or screw attachment when ordering.
- **Multi-ribbon option:** suitable for mounting on two-or-more slot metal ribbon holders.
- **Recognizable symbol:** clear identification for volunteer nurses of the CRI.

How to use / When it's useful

Quick answer: pin or mount the ribbon on your uniform or on a metal ribbon holder using the chosen clamp or screw attachment. In most cases, volunteers wear this ribbon on the uniform breast or group multiple ribbons on a single base for formal display.

Use the flat version for a standard, streamlined look or the Elite domed variant for a slightly more pronounced finish. Choose clamp attachment for easy removal or screw attachment for a more secure fit when mounting to a metal base or display.

Care and notes

Simple care keeps the ribbon presentable and suitable for daily use or display.

- Store in a dry place to avoid tarnish on metal parts.
- When ordering, confirm clamp or screw attachment to match your uniform or display needs.
- Only available in the 37 mm regulation format; other sizes are not offered.
- Limited availability: stock may run low, so check purchase options before finalizing.

FAQ

Q: What size is the **nastrino Infermiere Volontarie CRI**?

A: The ribbon is available only in the regulation 37 mm size. This is the standard format listed for this product.

Q: Which finishes are available?

A: Two finishes are offered: a flat regulation finish and a slightly domed "Elite" variant for a



different visual profile.

Q: Can I choose how the ribbon attaches to my uniform?

A: Yes. You can select a rear clamp attachment for easy removal or a screw attachment for a more permanent fit when ordering from the dropdown.

Q: Can multiple ribbons be mounted together?

A: Yes. For mounting two or more ribbons on a single metal base, choose an appropriate multi-slot ribbon holder designed to house multiple bars.

Q: Is this ribbon suitable for tactical or duty uniforms?

A: The ribbon is designed for CRI volunteer nurses and can be used on uniforms or display holders; it is compatible with duty and uniform displays commonly found in tactical or military surplus collections.

Q: Is stock limited?

A: Yes, availability may be limited. If you need the ribbon for an upcoming event or as part of a collection, we recommend ordering promptly.

Tip: Choose the finish and attachment that best fit your uniform or display—flat for a traditional look, Elite for a subtle domed effect; select clamp for easy removal or screws for a secure mount, and browse related ribbon holders to complete your setup.

Product attributes:

Nastrini : Formato 37 mm., Confezione Regolare, Confezione Elite, Confezione Regolare con Attacco a Clip, Confezione Regolare con Attacco a Vite, Confezione Elite con Attacco a Clip, Confezione Elite con Attacco a Vite, 50 cm. Nastro 37 mm.

Product gallery:

