



Dental Officer Service Ribbon



Product price:

€8.00 tax included

Product codes:

Reference: -

EAN13: 0658606555945

UPC: 658606555945

Product features:

Product short description:

The "Nastrino Ufficiale Odontoiatra" is an elegant, distinctive ribbon accessory designed for dental professionals. Crafted from high-quality materials, it adds a refined touch to uniforms and work attire, suitable for conferences, professional events, or daily clinical use.

- **Material:** Made from durable, high-quality fabric for a neat, professional finish.
- **Use:** Suitable for conferences, ceremonies, events and everyday wear in the dental office.
- **Design:** Refined, discreet style that complements uniforms without being obtrusive.
- **Benefit:** Displays professional affiliation with subtlety and enhances your presentation.

Product description:

Product overview

The Nastrino Ufficiale Odontoiatra is a prestige ribbon intended for personnel who hold the title Ufficiale Odontoiatra. It is a small, wearable insignia designed for display on uniforms, presentation boards, or multi-ribbon holders.

Definition: a nastrino is a ribbon-style duty insignia that identifies official rank or recognition. Suitable for military ribbons, duty insignia collections, and official ceremonial use within dental corps.



Main features

Clear, durable construction and traditional design make this nastrino a recognizable mark of rank and service.

- **Size:** 10 mm height × 37 mm width.
- **Design:** white background with a 2 mm orange horizontal stripe and a central round white insert bearing a red cross.
- **Materials:** made from high-quality, durable materials for long-lasting use.
- **Packaging options:** available in the standard flat format or the slightly domed "Elite" variant.
- **Rear attachment:** choose between clamps or screws for secure fastening.
- **Compatibility:** can be mounted on a multi-ribbon metal base (porta nastrini) to group multiple ribbons.

How to use / When it's useful

Quick answer: wear the nastrino on the uniform or mount it on a metal ribbon bar to show the Ufficiale Odontoiatra title at ceremonies and official events. It provides a compact, formal way to display professional recognition.

In most cases, choose the flat format for standard uniform use and the "Elite" domed variant for presentations or display boards. Select the rear attachment (clamps or screws) based on whether you need a removable or more permanent fit.

Care and notes

Basic care extends the nastrino's appearance and function; handle attachments carefully and store ribbons flat or on a holder when not in use.

- Clean gently with a soft, dry cloth; avoid solvents or abrasive cleaners.
- Check clamps or screws periodically to ensure a secure fit.
- Choose the correct rear attachment at purchase—clamps for easy removal, screws for a firmer mount.
- Limited availability: this item is in high demand, so allow lead time when ordering.

FAQ

Q: What are the dimensions of the Nastrino Ufficiale Odontoiatra?

A: The ribbon measures 10 mm in height and 37 mm in width. These are the standard dimensions for this nastrino.

Q: What rear attachment options are available?



A: You can choose between clamps or screws. Both options provide stable fastening; clamps allow easier removal while screws give a more permanent fit.

Q: Is the nastrino available in different formats?

A: Yes. It is offered in the standard flat format and a slightly domed "Elite" variant for a different visual effect.

Q: Can I display multiple nastrini together?

A: Yes. Multiple ribbons can be mounted on a metal multi-ribbon base (porta nastrini) to create an organized display.

Q: How can I personalize my nastrino?

A: Personalization options are limited to selecting the format (flat or Elite) and the rear attachment type (clamps or screws) at the time of purchase.

Q: What materials is the nastrino made from?

A: The nastrino is made from high-quality, durable materials intended for long-term wear and display. Specific material details are provided by the supplier if required.

Tip: When ordering, select the preferred rear attachment and format to match your uniform or display needs; check related ribbon holders and accessories to complete your presentation without rushing the purchase.

Product attributes:

Nastrini : Formato 37 mm., Confezione Regolare, Confezione Elite, Confezione Regolare con Attacco a Clip, Confezione Regolare con Attacco a Vite, Confezione Elite con Attacco a Clip, Confezione Elite con Attacco a Vite

Product gallery:

